Phone: +91-9896600003

Email: sales@naugralabequipments.com

#### **Product Name:**

Pharma Automatic Tablet and Capsule Dissolution Tester

### **Product Code:**

NLAB-PHARMACYLABM11001



## **Description:**

Pharma Automatic Tablet and Capsule Dissolution Tester

## **Technical Specification:**

It has four vessels and three poles ,arranged on one line ;

The head part of the instrument can be artificially turned over smoothly and flexibly .

The basket units and paddle units are made of imported stainless steel (ss316L)

The magnetism-pumped circular water current can heat the system evenly ,and the bath liquid can achieve an equal temperature.

A MPU is used in dissolution tester to automatically control the temperature, rotational speed and time, and the parameters can be preset ant any time.

Temperature range: ambient to +39.9°C

Speed range: (20-200)rpm

Speed accuracy: ±2rpm

Temperature Stability: ±0.3°C

Phone: +91-9896600003

Email: sales@naugralabequipments.com

Automatic time range: (5-900)min

Time accuracy: ±0.5min

Paddle radial run-out: ±0.5min

Basket radial run-ou: ±0.1min

Power: 220v/50HZ/1000W or 110v/60HZ/1000W

# Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com **Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003