

|                                                                                                                                                                                  |                                               |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
| <b>Product Name :</b><br>Plastic Needle Cutter                                                                                                                                   | <b>Product Code :</b><br>NLABMEDICAL-M3370001 |
|                                                                                                                                                                                  |                                               |
| <b>Description :</b>                                                                                                                                                             |                                               |
| Plastic Needle Cutter                                                                                                                                                            |                                               |
| <b>Technical Specification :</b>                                                                                                                                                 |                                               |
| To Terminate Needles, insert used needle fully into small insert used needle fully into small orifice and hold squeeze handles with a fast squeezing action until needle is cut. |                                               |
| To Terminate Syringes, insert syringe nozzle fully into large orifice and hold squeeze handles with a fast squeezing action until nozzle is cut.                                 |                                               |

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003