

Product Name : Slide Mounter	Product Code : NLAB-EDUSCHOOLEQ29020
Description : Slide Mounter	
Technical Specification : A Slide Mounter can be used for several things. If you've got a diatom test slide from Klaus Kemp you will note a small "finder ring" around the 8 diatoms to help locate them quickly. In addition, the edge of the round coverslip is sealed to the slide very neatly, again using a ringer. On coverslips that I ring, I let the dark colored edge sealant dry, then add a thin white ring on top to allow identifying the top of the slide while not wearing glasses. Another use is making support rings to prevent crushing a specimen when applying the cover slip. This was my original motivation since I wanted to make whole mounts of deer ticks without mashing them in the process. I tried using punched paper (lots of stray	

fibers) and punched tape (melted by xylene in the mountant) with poor results;
shellac rings are easy to make and work very well.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003