

Product Name : Bile Salts Brilliant Green Starch Agar	Product Code : NLAB-CHEMICALSIND21035
Description :	
Bile Salts Brilliant Green Starch Agar	
Technical Specification :	
Bile Salts Brilliant Green Starch Agar	
Bile Salts Brilliant Green Starch Agar is recommended for the selective isolation and identification of Aeromonas hydrophila from food and environmental specimens.	
Directions	
Suspend 45 grams in 1000 ml distilled water. Heat to boiling with occasional agitation to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.	
Quality Control	
Appearance: Light yellow to greenish yellow homogeneous free flowing powder	
Reaction: Reaction of 4.5% w/v aqueous solution at 25°C. pH : 7.2±0.2	
Gelling: Firm, comparable with 1.5% Agar gel	
pH: 7.00-7.40	
Storage: Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C.	
Colour and Clarity of prepared medium: Green coloured, slightly opalescent gel forms in Petri plates.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003